



ResMed

Actual Size: 5 1/2" (Width)
x 6 7/8" (Height)

Ruler is for scaling purposes.

How to use the fitting template:

1. Please ensure the template is printed to scale
2. Cut out along dotted line
3. Refer to diagram on template to determine appropriate cushion size

AirFit™ N20

Nasal mask
Nasal mask for Her



3.7 - 4.5 CM
1 1/2 - 1 3/4 IN



< 3.7 CM
< 1 1/2 IN



4.5 - 5.5 CM
1 3/4 - 2 1/8 IN



Rx Only 0123